

EXHIBIT K

Abstract
#12044***A priori* filtering of post-operative circulating tumor DNA predicts recurrence in post-metastasectomy colorectal cancer patients without knowledge of tumor genotype**Michael J. Overman¹, Ariel Jaimovich², Drew Kennedy², Michelle Tan², Danielle Gavino², Stefanie Mortimer², Darya Chudova², AmirAli Talasaz², Justin Odegaard², Scott Kopetz¹¹The University of Texas MD Anderson Cancer Center, Houston, TX, USA; ²Guardant Health, Inc., Redwood City, CA, USA

THE UNIVERSITY OF TEXAS

MD Anderson
Cancer Center**Introduction**

- ctDNA in post-operative colorectal cancer (CRC) patients correlates with molecular residual disease and may be useful for prognostication and to guide adjuvant therapy decision making.^{1,2}
- We previously demonstrated that post-operative ctDNA is strongly associated with disease recurrence in patients with metastatic CRC undergoing curative intent surgery ($p=0.004$).³
- Initial studies employed clinically impractical assays indexed to individual patient specific tumor tissue-derived mutations or were confounded by non-tumor-associated somatic alterations, including variants related to clonal hematopoiesis¹.
- We previously demonstrated that using a highly sensitive CRC next-generation sequencing (NGS) panel, the detection of post-operative ctDNA does not require foreknowledge of known somatic alterations.³
- We developed a variant classifier to expand on this ctDNA only approach and to further differentiate tumor-derived alterations from non-tumor derived alterations with the goal of increasing specificity of ctDNA detection in post-operative CRC patients.

Methods



- CRC patients planned for hepatic metastasectomy were prospectively enrolled in an IRB approved trial (LAB10-0982).
- Pre-operative and post-operative plasma was sequenced to high depth using a 23-gene NGS panel with 96% theoretical sensitivity for CRC.³
- 51 metastatic colorectal cancer patients with both pre and post ctDNA results were recruited at a single institution (Tables 1,2). Tumor tissue was sequenced using this panel or local testing.
- ctDNA profiles from 4000 CRC pts (Guardant Health, Redwood City, CA) were used to train a variant classifier to exclude non-tumor derived alterations.
- The variant classifier was designed to identify ctDNA mutations that originate from the tumor, differentiating them from non-tumor derived mutations.

Table 1: Cohort demographics

Number of unique patients	51
Median age at diagnosis (range)	55 years (33-76)
Gender	60.8% Male 39.2% Female
Histological Grade	98% Moderately Differentiated 2% Poorly Differentiated
Primary site	21.6% Right-sided 78.4% Left-sided
Presentation	15.7% Metachronous 84.3% Synchronous

Table 2: Cohort clinical features

Neoadjuvant chemotherapy	80.4%
Median number of resected tumors	2
Lymph node positive primary	66.7%
KRAS mutation positive	43%
Median time surgery to post-operative sample (range)	18 days (13-123 days)
Median follow-up (range)	42.7 months (4.4 – 59.4 months)
Recurrence	72.5% (37 patients)
Median Time to Recurrence (range)	7.8 months (1.2 – 34.5 months)

Results**Variant Allele Frequency distribution****Conclusions**

- Recurrence prediction using post-operative somatic variant detection alone is fraught by a high clinical false positive rate.
- Many non-tumor derived mutations occur at low variant allele frequencies. However, relying solely on an allele frequency threshold to differentiate between tumor derived and non-tumor derived mutations would exclude many clinically relevant mutations.
- Filtering using tumor tissue is effective but may be clinically impractical due to added complexity and cost.
- Filtering using a novel variant classifier, without foreknowledge of tumor genotype eliminated false positives while maintaining clinically acceptable sensitivity.
- *A priori* variant classification may enable clinically feasible ctDNA diagnostics for adjuvant decision making in early-stage disease.

References

1. Tie J, et al. (2016). Circulating tumor DNA analysis detects minimal residual disease and predicts recurrence in patients with stage II colon cancer. *Sci Transl Med* 8(346)
2. Diehl F, et al. (2008). Circulating mutant DNA to assess tumor dynamics. *Nat Med* 14(9)
3. Overman, et al. (2017). Circulating tumor DNA (ctDNA) utilizing a high-sensitivity panel to detect minimal residual disease post liver resection and predict disease recurrence. *JCO* 35(suppl).